



## ■ 最大積載表

### -Universal Trays- 最大積載数

	Size	KM Mini 200 x 295mm	KM 356 x 300mm	KS 430 x 300mm	SM 560 x 400mm	VKS 760 x 600mm
Spring clamps	10ml	23	58	68	117	280
	25ml	15	29	34	64	138
	50ml	15	27	27	63	136
	100ml	15	24	24	49	114
	250ml	6	12	15	23	52
	500ml	4	8	9	15	35
	1000ml	3	5	5	11	21
	2000ml	1	2	2	6	12
	3000ml	-	-	1	4	9
	5000ml	-	-	1	2	6
Test tube racks (with hinged foot)		1	3	4	5	8

### -Rack Systems with Clamping Strips- 最大積載数

	Size	KM Mini 200 x 295mm	KM 300 x 300mm	KS 400 x 275mm	SM 560 x 400mm	VKS 760 x 600mm
Erlenmeyer flasks	10ml	20	-	-	-	-
	25ml	16	-	-	-	-
	50ml	12	20	25	49	104
	100ml	9	12	20	42	82
	250ml	6	9	12	20	45
	500ml	4	4	6	12	32
	1000ml	2	2	4	6	20
	2000ml	1	1	2	6	12
	3000ml	-	1	1	4	9
	5000ml	-	1	1	2	6
Test tube racks		1	1	2	5	8
Separating funnels	100ml	-	4	4	10	16
	250ml	-	3	3	6	12
	500ml	-	-	2	4	8
	1000ml	-	-	1	4	6
	2000ml	-	-	1	2	4

